



PRACTICAL INNOVATION

Ziath

# ZIATH™ MANUAL HEAT SEALER

Product ID – ZMHS-01



The **Ziath Manual Heat Sealer** has been designed to be compact, easy to use and ergonomic. Sealing a wide range of plates via heat sealing is quick and simple. The sealing protects precious samples from evaporation and contamination when performing PCR or when storing samples. The device runs on conventional power supply, and therefore needs no compressed air.

[www.ziath.com](http://www.ziath.com)


## System Specification

---

- Ideal entry-level sealer
- Benchttop manually operated plate sealing device
- Heat-seals SBS format micro well plates
- Adjustable plate heights
- Seal temperature  $170^{\circ}\text{C} \pm 2^{\circ}\text{C}$
- Warm up time less than 10 minutes
- Ergonomic pull-down operation
- LED displays for system on/system ready for use
- Lightweight system — only 6kgs (14lbs)
- 110/220 vac power compatible
- No compressed air required
- Consistent, single action sealing
- Heat sealing foils available from Ziath

## System Footprint

---



**For more information please contact Ziath Ltd.**

[www.ziath.com](http://www.ziath.com) • [sales@ziath.com](mailto:sales@ziath.com)

**Tel +44 (0)1223 655577 • Fax +44 (0)1223 790069**