



PRACTICAL INNOVATION

Ziath

# ZIATH™ OPTILABELER

Product ID – ZOL-01



The **Ziath OptiLabeler** makes light work of labeling vials and tubes. This device is a simple heavy-duty semi-automatic electric label dispenser for labels from 6mm to 50mm widths to any length. The dispenser advances at over 10cm per second and allows you to peel small labels with ease.

[www.ziath.com](http://www.ziath.com)


## Device Features

---

- Ideal for applying labels on tubes and vials
- Will dispense labels and die-cut parts
- Will dispense labels from 6mm to 50mm width
- Will dispense labels from 6mm to 75mm length
- Will accommodate core sizes from 25mm up to 75mm
- Photosensor accuracy
- Automatically removes paper, vinyl, acetate, polyester, foil, etc. from its liner
- Automatically rewinds liner scrap for a cleaner working environment
- Liner scrap is easily removed
- No tools required for adjusting to different size labels
- All metal construction
- No maintenance required

## Device Footprint

---



**For more information please contact Ziath Ltd.**

[www.ziath.com](http://www.ziath.com) • [sales@ziath.com](mailto:sales@ziath.com)

**Tel +44 (0)1223 655577 • Fax +44 (0)1223 790069**