

If you would no longer like to receive emails from Ziath then please click [here](#).
If you cannot see the email, please view it on our website [here](#).

Ziath


Practical Innovation

The Ziath Newsletter

www.ziath.com
10th October 2009

New Heat Sealer

[Click to contact us](#)



Ziath have teamed up with the leading manufacturer of heat sealers and are pleased to launch a new compact semi-automated heat sealer.

- Seal SBS format micro well plates
- Seal plate height with no adjustment from 9 - 48mm high
- Temperature range from ambient to 200 degrees Celsius
- Seal time from 0.1 - 9.9 seconds
- Simple slide operation
- Combined temperature and seal time display
- Light weight system only 6kgs (14lbs)
- 110/220 vac power compatible
- No compressed air required